

MICRO APPLICATORS

DISPOSABLE MICRO APPLICATORS



CONFIGURATION OVERVIEW

- Material:** Medical level Polypropylene handle, Nylon Brush
- Type:** Non-Sterile
- Description:**
- Ideal for application of sealants, bonding agents, etchants, hemostatic solutions, cavity liners and disclosing solutions.
 - Non-absorbent, non-linting fibers, suspend solutions as small as 1/8 of a drop without spilling, dripping or waste.
 - The bendable neck stays in place to facilitate the placement of material in difficult areas.
 - Flexible plastic handle to pinpoint accuracy application of any liquids.





Class 1 Medical Device according to EU MDR 20017/745 according to rule 5 Annex VIII






COMPLIANT WITH: EN ISO 14971:2019, EN ISO 15223-1:2021, EN ISO 10993-1:2018, EN ISO 10993-5:2009, EN ISO 10993-10:2010, EN ISO 20417:2021


THECNICAL SPECIFICATION

Size	Regular	Fine	Superfine
Length of the applicator (including the fiber tip) (mm) ± 0.1 mm	102	101.5	100.2
Diameter of the fiber tip (mm) ± 0.3 mm	3	1.9	1.5
Length of the fiber tip (mm) ± 0.5 mm	3.5	3.5	2.5

PRODUCT VARIANTS

Regular	Art-No.	MICAP25-AKZ-BL-REG	MICAP25-AKZ-LI-REG	MICAP25-AKZ-FG-REG	MICAP25-AKZ-PI-REG
	Color	Med Blue	Lilac	Fresh Green	Pink
					

Fine	Art-No.	MICAP25-AKZ.YE.FNE	MICAP25-AKZ-MI-FNE	MICAP25-AKZ-MC-FNE	MICAP25-AKZ-LU-FNE	MICAP25-AKZ-CE-FNE
	Color	Yellow	Mint Green	Magic Blue	Laguna Blue	Crispy cedro
						

Superfine	Art-No.	MICAP25-AKZ-WH-SFN
	Color	White
		

SHELF LIFE : 5 years

PACKAGING

Pcs/dispenser	100
Dispenser/box	4
Boxes/carton	40
Carton size	52x42x24.5 cm
Carton weight	16 Kg

STORAGE AND HANDLING

In dry, dark, dust-free and ventilated space at room temperature. Keep away from direct sunlight.

INSTRUCTIONS FOR USE

Choose the suitable size of applicators to be coated with medicine to be applied to the affected area.

WARNING

- The fiber brush cannot be used for a long time with absorbed large volume of chemical solution. Avoid long usage to prevent the fiber brush to fall out.
- Do not use with reagents containing BUTANONE.



Rev. 3 23-12-2024