

SUNFLOWER LATEX GLOVES

SINGLE USE DISPOSABLE LATEX EXAMINATION GLOVES

CONFIGURATION OVERVIEW

Features

- Cuff strengthened by a rolled cuff
- Easy donning and removal
- Excellent tactile sensibility (Micro textured fingertips)
- Suitable for food handling purposes
- Ambidextrous

Material: Natural rubber latex

Type: Non-Sterile
Powder-free
Polymer coated



TECHNICAL SPECIFICATIONS

Physical Properties:

Force at Break (%)before and after aging	≥ 6.0 N
Freedom from holes	AQL ≤ 1,5
Residual Powder	Max. 2.0 mg/glove
Protein Content	Max. 50 µg/g
Tensile Strength	Before aging: Min 18 Mpa After Aging: Min. 14 Mpa
% Elongation	Before Aging: Min. 650% After Aging: Min 500%

Location of Thickness Measurement	Single Wall (mm)
Finger Tip	≥ 0.08
Palm	≥ 0.08

Available Length: Med. ≥ 240 mm
Available Sizes /Palm width (mm) : XS : ≤80 - S: 80±10 - M: 95±10 - L: ≥110±10
Weight M size: 5.9 g ± 0.3 g
Shelf life: 5 years

QUALITY CERTIFICATION: Class I Medical Device as per MDR 2017/745
Eu Reg. 2016/425 PPE CAT III

Conform to: EN 455-1:2020+A1:2022, EN 455-2:2015, EN 455-3:2023, EN 455-4:2009, EN ISO 15223- 1:2021, EN ISO 21420:2020, EN ISO 20417:2021, EN ISO 14971:2019/A11:2021, EN ISO 10993-1:2018, EN ISO 10993-5:2009, EN ISO 10993-10:2010, EN 374-1:2024, EN 374-2:2019, EN 374-4:2019, EN 374-5:2016, EN 16523-1:2015+A1:2018, ASTM F 1671:2013, ASTM D3578:2019, ASTM D7161:2016, ASTM D7160:2016, EC Reg. 1935/2004 and its subsequent amendments

PRODUCT VARIANTS

Art-No.	GLODL29-AKZ-WH-XS1/S01/M01/LL1
Color	White

Available Sizes: XS S M L

Tel.: +41 (0) 91 921 14 92
 Web: akzenta.com
 Mail: business@akzenta.com

Visit us at
akzenta.com

USt. ID Nr. DE815221533
 VAT No. CHE-337.693.284 IVA
 Company No CHE-337.693.284

PACKAGING

Pcs/Dispenser:	100
Dispenser/Carton:	10
Carton Size:	24.9 X 34 X 25 cm
Carton weight:	XS: 5,74 Kg – S: 6.09 kg – M: 6,82 Kg - L: 7.5 Kg
Carton / layer:	9
Layer / pallet:	7
Carton / pallet:	63
Pallet size:	80x120x186cm
Pallet weight :	XS: 382 Kg. - S: 404 Kg. - M: 449 Kg. - L: 493 Kg

PRODUCT INSTRUCTIONS FOR USE

- 1) Before usage, inspect the gloves for any defects or imperfections.
- 2) **Usage is for single use only.** If gloves are reused, the risk of contamination and infection increases due to improper cleaning processes, as well as an increased risk of holes and tears during reuse caused by weakening of the gloves.
- 3) Sizing: Select the right size of glove for your hands.
- 4) Donning: Hold the glove by the cuff with one hand. Align the glove thumb with your other hand's thumb and slide your hand into the glove, inserting one finger into each glove finger. Pull on the glove palm to achieve a good fit. Don the other glove using the same procedure.
- 5) Inspection: punctures or tears may occur after donning. Inspect the glove after donning and stop using them immediately if any damage is found.
- 6) Doffing: Hold the glove cuff and pull towards the fingers until the glove comes off.
- 8) Disposal: For the proper disposal of all used gloves, please follow your institution's policies for disposal.

STORAGE: Store in cool, dry place at room temperature. Opened box should be shielded from exposure to direct sun or fluorescent/ultraviolet light sources. Storage above 30° will lead to accelerated aging and should be avoided.

