

PRIME|ONE OAT COATED NITRILE GLOVES

SINGLE-USE DISPOSABLE NITRILE EXAMINATION GLOVES

CONFIGURATION OVERVIEW

| | |
|------------------|--|
| Features | Cuff strengthened by a rolled cuff Easy donning and removal Excellent tactile sensibility (Micro textured fingertips) Coated with colloidal Oatmeal |
| Shape: | Ambidextrous |
| Material: | Nitrile |
| Type: | Non-Sterile Powder-free Latex-free Chlorinated |



TECHNICAL SPECIFICATIONS

Physical Properties:

| | |
|---------------------------------------|------------------|
| Elongation at Break (%) | ≥500% |
| Force Break (N) Before & After ageing | ≥6N |
| Colloidal Oatmeal content | ≥ 5 mg/glove |
| Location of Thickness Measurement | Single Wall (mm) |
| Fingertip | ≥ 0.09 ± 0.02 mm |
| Palm | ≥ 0.07 ± 0.02 mm |

| | |
|--------------------------------|----------------------------|
| Available Length: | minimum 240 mm |
| Median Palm Width ± 4mm | XS : 76 S: 86 M: 98 L: 107 |
| Weight M size: | 3.2 g ± 0.2 g |

SHELF LIFE: 5 years

QUALITY ASSURANCE :


Class I Medical Device under MDR EU REG. 2017/745
Category III Personal Protective Equipment EU Reg. 2016/425

CONFORM TO: EN ISO 374-1:2016+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 16523-1:2015+A1:2018, EN ISO 21420:2020, ISO 16604:2004, EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, EN 455-4: 2009, EN ISO 14971:2019/A11:2021, EN ISO 10993-1:2021, EN ISO 10993-5:2009, EN ISO 10993-10:2010, EN ISO 10993-11:2018, ISO 2859-1:1999/Amd1:2011, ISO 15223-1:2021/Amd.1:2025, ASTM D6319, ASTM D7160

Conform to : EC REG. 1935/2004



PRODUCT VARIANTS

| | |
|---------|---|
| Art-No. | GLOBN53-AKZ-DW-XS1/S01/M01/LL1 |
| Color | Dawn Blue  |

Available Sizes:

XS S M L

PACKAGING:
Pcs/Unit: 100

Unit/Carton: 10

Carton/Pallet: 84

Carton/Layer: 12

Layer/Pallet: 7

Carton Dimension: 26.2x25.2x25.7 cm

Carton weight : XS 3,4 Kg. – S 3.7 Kg. – M 4.0Kg. – L 4.3 Kg.

Pallet Dimension: 80x120x200

Pallet/Weight: XS 305 Kg. – S 330Kg. – M 356 Kg. – L 382 kg.

PRODUCT INSTRUCTIONS FOR USE

- 1) Before usage, inspect the gloves for any defects or imperfections.
- 2) **Usage is for single use only.** If gloves are reused, the risk of contamination and infection increases due to improper cleaning processes, as well as an increased risk of holes and tears during reuse caused by weakening of the gloves.
- 3) Sizing: Select the right size of glove for your hands.
- 4) Donning: Hold the glove by the cuff with one hand. Align the glove thumb with your other hand's thumb and slide your hand into the glove, inserting one finger into each glove finger. Pull on the glove palm to achieve a good fit. Don the other glove using the same procedure.
- 5) Inspection: punctures or tears may occur after donning. Inspect the glove after donning and stop using them immediately if any damage is found.
- 6) Doffing: Hold the glove cuff and pull towards the fingers until the glove comes off.
- 8) Disposal: For the proper disposal of all used gloves, please follow your institution's policies for disposal.

STORAGE: Store in cool, dry place at room temperature. Opened box should be shielded from exposure to direct sun or fluorescent/ultraviolet light sources. Storage above 30° will lead to accelerated aging and should be avoided.

